

TruePure™ dNTPs

About us

Since 2001, Canvax has been an original Manufacturer & Supplier of the Most Innovative, **Enzymes, Kits, Buffers and Reagents**. Due our customizable and adaptable manufacturing structure, reduced overhead and capacity to BULK manufacturing with the highest quality standards, we are **leading the Life Science industry in terms of product performance, manufacturing capabilities and customer-service**.

Quality & Certifications

As a leading reagent manufacturer, we are proud of our **first-class manufacturing facilities**, located in Valladolid (Spain) with a total area of 3,400m², specifically **designed for production Enzymes, Buffers and Kits on commercial scale**. Canvax operates in its first-class facilities under certified:

- >> **ISO 9001:2015**.
- >> **ISO 13485:2016** for the Research use Raw Materials and Ancillary Reagents.
- >> **cGMP**, current Good Manufacturing Practice and in continuous improvement according to ICH Q7.

TruePure™ dNTPs

The need of ultrapure deoxynucleotide triphosphates (dNTPs) in long PCRs, RT-PCRs and proofreading amplifications is critical. In this line, Canvax™ offers **TruePure™ dNTPs**, manufactured by a Certified Quality Management System where each dNTP lot is assayed by the most stringent criteria regarding Purity and Functionality.

TruePure™ dNTPs has the following **Advantages & Features**:

- >> **Highest Purity**: greater than 99% confirmed by HPLC.
- >> **Highly stable**.
- >> **Proven Performance for PCR** and most common Molecular Biology applications.
- >> **Time-saving**: ready-to-use format that eliminates experiment preparation time.
- >> **High Quality**: free of DNase, Protease, Phosphatase, Nuclease, Human or E. coli DNA contamination.



Ordering info:

Description	Concentration of each nt	Total Concentration	Cat. No.	Size
TruePure™ dNTP Mix	2 mM	8 mM	N0030-S	1 mL
	2.5 mM	10 mM	N0031-S	1 mL
	25 mM	100 mM	N0032	1 mL
	10 mM	40 mM	N0033-S	1 mL
TruePure™ dNTP Set		100 mM	N0098-L	4 x 1 mL

Canvax Reagents S.L.U. | A: c/Luis de Mercado 19, Parque Tecnológico de Boecillo, PC: 47151, Boecillo; Valladolid, Spain
T: +34 983 54 85 63 | M: info@canvaxbiotech.com | W: canvaxbiotech.com

Prices valid until December 31th, 2023. Taxes and transport not included.

For Research Use Only (RUO Mark). Not for use in Diagnostic procedures.

© 2023 Canvax Reagents, S.L.U. All Rights Reserved. The **tradenames** of Canvax products referenced herein, whether or not appearing in bold characters or with the trademark symbol ® or ™, are trademarks of entities of Canvax. The names of actual companies and products mentioned herein maybe the trademarks of their respective owners.