

Oxford Accupet EVO

Single Channel Pipette

Volume lock setting to prevent accidental volume drifts

Made using highly durable material for corrosion and physical shock resistance

Highly durable PVDF tip cones


Corrosion resistant tip ejector with a unique shock absorbing mechanism



Magnetically assisted piston for precise and consistent results



Innovative spring and seal design for one of the lowest plunger forces for a fatigue free experience



Universal tip compatibility

Model#	Description	Volume	No. of Channels
AE-2.5	Accupet EVO Single Channel Adjustable Vol. Pipette	0.1-2.5 μ l	Single
AE-10	Accupet EVO Single Channel Adjustable Vol. Pipette	0.5-10 μ l	Single
AE-20	Accupet EVO Single Channel Adjustable Vol. Pipette	2-20 μ l	Single
AE-50	Accupet EVO Single Channel Adjustable Vol. Pipette	5-50 μ l	Single
AE-100	Accupet EVO Single Channel Adjustable Vol. Pipette	10-100 μ l	Single
AE-200	Accupet EVO Single Channel Adjustable Vol. Pipette	20-200 μ l	Single
AE-1000	Accupet EVO Single Channel Adjustable Vol. Pipette	100-1000 μ l	Single
AE-5000	Accupet EVO Single Channel Adjustable Vol. Pipette	500-5000 μ l	Single
AE-10000	Accupet EVO Single Channel Adjustable Vol. Pipette	1000-10000 μ l	Single

3 Year Warranty

For more details, visit us at www.oxfordlp.com

AN ISO 9001 : 2015 & ISO 13485 : 2016 CERTIFIED MFG. COMPANY

*OEM request can be entertained. Specifications can be changed without notice for quality improvement.

Oxford Accupet EVO

Multichannel Pipette

Volume lock setting to prevent accidental volume drifts

Fully autoclavable

Ergonomic design with proper weight distribution



Highly durable PVDF tip cones



Spring-loaded tip cone for less strain & proper tip fitment for precise pipetting



Color coding for easy identification



Large and clear 4-digit display

Model#	Description	Volume	No. of Channels
AE8-10	Accupet EVO 8-Channel Adjustable Vol. Pipette	0.5-10 μ l	8
AE8-20	Accupet EVO 8-Channel Adjustable Vol. Pipette	2-20 μ l	8
AE8-50	Accupet EVO 8-Channel Adjustable Vol. Pipette	5-50 μ l	8
AE8-100	Accupet EVO 8-Channel Adjustable Vol. Pipette	10-100 μ l	8
AE8-200	Accupet EVO 8-Channel Adjustable Vol. Pipette	20-200 μ l	8
AE8-300	Accupet EVO 8-Channel Adjustable Vol. Pipette	30-200 μ l	8
AE12-10	Accupet EVO 12-Channel Adjustable Vol. Pipette	0.5-10 μ l	12
AE12-20	Accupet EVO 12-Channel Adjustable Vol. Pipette	2-20 μ l	12
AE12-50	Accupet EVO 12-Channel Adjustable Vol. Pipette	5-50 μ l	12
AE12-100	Accupet EVO 12-Channel Adjustable Vol. Pipette	10-100 μ l	12
AE12-200	Accupet EVO 12-Channel Adjustable Vol. Pipette	20-200 μ l	12
AE12-300	Accupet EVO 12-Channel Adjustable Vol. Pipette	30-300 μ l	12

3 Year Warranty

For more details, visit us at www.oxfordlp.com

AN ISO 9001 : 2015 & ISO 13485 : 2016 CERTIFIED MFG. COMPANY

*OEM request can be entertained. Specifications can be changed without notice for quality improvement.